ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 800-927-4669
Tel: +1 903-675-9321
www.argonmedical.com

EMERGO EUROPE

Westervoortsedijk 60
6827 AT Arnhem
+31 70 345 8570

Emergo Consulting (UK) Limited
c/o Cr 360 - UL International
Compass House, Vision Park Histon
Cambridge CB24 9BZ
United Kingdom

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CH REP
MedEnvoy Switzerland
Gotthardstrasse 28
6302 Zug
Switzerland

<https://www.argonmedical.com/resources/product-information>

SSCP Available: <https://www.argonmedical.com/resources/product-information>

For a copy of this device's current European Summary of Safety and Clinical Performance (SSCP), please go to the European database on medical devices (Eudamed), where it is linked to the basic UDI-DI. <https://ec.europa.eu/tools/eudamed>.

| | |
|----------------------|--|
| English | The symbols glossary is located electronically at www.argonmedical.com/symbols |
| Brazilian Portuguese | O glossário dos símbolos está localizado eletronicamente em www.argonmedical.com/symbols |
| Bulgarian | Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols |
| Czech | Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols |
| Danish | Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols |
| German | Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden |
| Greek | Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols |
| Spanish | El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols |
| Estonian | Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols |
| Finnish | Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols . |
| French | Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols |
| Hungarian | A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols |
| Italian | Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols |
| Lithuanian | Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols |
| Norwegian | Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols |
| Dutch | Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols |

| | |
|------------|--|
| Polish | Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols |
| Portuguese | O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols |
| Romanian | Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols |
| Slovak | Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols |
| Swedish | Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols |
| Turkish | Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır |
| Arabic | يتوفر مصدر الرموز إلكترونياً على www.argonmedical.com/symbols |

IFU1100M
Date of Issue: 11/21/2024 Rev A

en - English

WORKER™ Guidewire (Standard, Amplatz and Bentson)

Intended Use/Purpose

Guidewires are intended to facilitate the percutaneous placement of intravascular and non-vascular devices during diagnostic and intervention procedures.

Device Description

WORKER™ guidewires are made of stainless steel core with a polytetrafluoroethylene (PTFE) coated stainless steel spring wire that helps reduce friction during procedure. The guidewires are compatible with medical devices with a lumen larger than the diameter size listed on the label and shorter than the labelled length. Physician judgement required to select the appropriate guidewire compatible with other medical devices to be used.

Indications for use

Guidewires are indicated for use in angiographic procedure to introduce and position catheters and interventional devices within the coronary and peripheral vasculature. Guidewires are also intended to facilitate the percutaneous placement of peripheral intravascular and non-vascular devices during diagnostic and intervention procedures. The Argon guidewires present an indirect benefit to the patient by enabling diagnosis or procedures targeting the heart and the central circulatory system, the peripheral vascular, or non-vascular procedures related to gallbladder or biliary obstruction and percutaneous drainage. The selection of guidewire is based on physician judgement based on the type of procedure being performed.

Duration

Transient, less than 60 minutes.

Intended Users

The product is intended for use by medical professionals and must only be used by qualified personnel who are familiar with the technique.

Contraindications for Use

There are no known contraindications.

Warnings

- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appears open or if the expiry date has been exceeded.
- Do not advance the wire against resistance until the cause of the resistance has been determined by fluoroscopy. Excess force against resistance may result in damage to guidewire, or catheter, or vessel perforation.
- Do not excessively torque the guidewire.
- Do not entrap or over-rotate the distal tip of the guidewire, it can lead to wire rupture.
- Do not excessively bend the guidewire, it may result in wire fracture.
- Guidewire placement should be monitored during placement and manipulation using fluoroscopic or suitable imaging method.
- Do not withdraw a guidewire through a needle. Straighten the guidewire in order to withdraw the needle.

Precaution

- Twisting or entanglement of the guidewires may occur when more than one wire is being used simultaneously. This can be prevented by carefully isolating and marking the proximal end of the wires.

Potential Complications/Adverse Events

Potential complications include, but are not limited to the following:

- Perforation of vessel or arterial wall
- Thrombus formation
- Puncture site hematoma
- Vasospasm
- Infection
- Vessel dissection

Preparation

- Make sure that the guidewire is undamaged.
- Loosen the straightener together with the guidewire.
- Inject saline with syringe in a vertical position.
- Make sure the product is compatible with other equipment.

Directions for Use

1. Insert the flexible end of the guidewire into the device. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device.
 - 1.1. If "J" wire is used, prepare "J" tip by sliding J-straightener over tip to straighten. Insert J-tip inside J-Straightener into device hub gently.
2. Advance the guidewire through the device.
3. Position the guidewire.
4. Hold the guidewire in place and withdraw the device.
5. Advance the catheter over the guidewire.

Storage

Store at controlled room temperature. Avoid extended exposure to light.

Disposal

After use, handle and dispose in accordance with end-user facility's policies and procedures concerning biohazardous materials and waste.

NOTE : In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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