**ARGON MEDICAL DEVICES, INC.**  
 1445 Flat Creek Road  
 Athens, Texas 75751 USA  
 Tel: +1 (903) 675 9321  
 Tel: +1 (800) 927 4669  
 www.argonmedical.com

**UK REP** Emergo Consulting (UK) Limited  
 c/o Cr 360 - UL International  
 Compass House, Vision Park Histon  
 Cambridge CB24 9BZ  
 United Kingdom

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**CH REP**  
 MedEnvoy Switzerland  
 Gotthardstrasse 28  
 6302 Zug  
 Switzerland

35-2000-10CM

Date of Issue: 11/21/2024 Rev A

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Bulgarian	Речникът на символите е наличен в електронен вид на адрес <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Czech	Slovník symbolů v elektronické formě je umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Danish	Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
German	Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Spanish	El glosario de símbolos puede consultarse en línea en la dirección <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Finnish	Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
French	Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Lithuanian	Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Latvian	Elektroniskā simbolu vārdnīca ir pieejama <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Norwegian	Symbolforklaringen befinner seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>

Portuguese	O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovenian	Slovar simbolov se nahaja v elektronski obliki na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
Turkish	Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır

## EN - English

### GUIDEWIRES

#### Intended Use/Purpose

Guidewires are intended to facilitate the percutaneous placement of peripheral intravascular and non-vascular devices during diagnostic and intervention procedures.

#### Device Description

Guidewires are designed to facilitate the passage of catheters and sheath introducers into the vasculature utilizing the Seldinger or Modified Seldinger techniques. The device is made with stainless steel and available with or without PTFE coating. The guidewires are compatible with medical devices with a lumen larger than the diameter size listed on the label and shorter than the labelled length. Physician judgement required to select the appropriate guidewire compatible with other medical devices to be used.

#### Indications for Use

Guidewires are indicated to facilitate the percutaneous placement of peripheral intravascular and non-vascular devices during diagnostic and intervention procedures. The Argon guidewires present an indirect benefit to the patient by enabling diagnosis or procedures targeting the peripheral vascular, or non-vascular procedures related to gallbladder or biliary obstruction and percutaneous drainage. The selection of guidewire is based on physician judgement based on the type of procedure being performed.

#### Duration:

Transient, less than 60 minutes.

#### Intended Users

The product is intended for use by medical professionals and must only be used by qualified personnel who are familiar with the technique.

#### Contraindications for Use

There are no known contraindications.

#### Possible Complications/Adverse Events

Potential complications include, but are not limited to the following:

- Perforation of vessel or arterial wall
- Thrombus formation
- Puncture site hematoma
- Vasospasm
- Infection
- Vessel dissection

#### Warnings

- This product is not intended for use in the coronary arteries.
- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package appears open or if the expiry date has been exceeded.
- Do not advance the wire against resistance until the cause of the resistance has been determined by fluoroscopy. Excess force against resistance may result in damage to guidewire, or catheter, or vessel perforation.
- Do not excessively torque the guidewire.
- Do not entrap or over-rotate the distal tip of the guidewire, it can lead to wire rupture.
- Do not excessively bend the guidewire, it may result in wire fracture.
- Guidewire placement should be monitored during placement and manipulation using fluoroscopic or suitable imaging method.
- Do not withdraw a guidewire through a needle. Straighten the guidewire in order to withdraw the needle.

#### Precautions

- Twisting or entanglement of the guidewires may occur when more than one wire is being used simultaneously. This can be prevented by carefully isolating and marking the proximal end of the wires.

#### CAUTIONS for Movable Core Guidewires

While the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

#### Preparation

1. Make sure the product is compatible with other equipment.

**Note:** If unable to remove guide wire from dispenser, separate straightener from hoop dispenser.

2. Flush the guidewire with sterile heparinized normal saline or a similar isotonic solution by connecting a syringe to the luer hub of the guidewire dispenser.

3. Inject saline with syringe in a vertical position.

#### Direction for Use

1. Insert the flexible end of the guidewire into the needle. Avoid overly vigorous manipulation of the guidewire to prevent damage of vessel or device. If "J" wire is used, prepare "J" tip by sliding J-Straightener over tip to straighten. Insert J-tip inside J-Straightener into the needle hub gently.

2. Advance the guidewire through the needle.

3. Position the guidewire.

4. Hold the guidewire in place and withdraw the needle.

5. Advance the intravascular device or catheter over the guidewire.

#### Movable Guidewires

To lengthen the flexible segment of a movable-core guidewire, gently withdraw the core. Do not advance movable-core within vasculature.

#### Disposal

After use, handle and dispose in accordance with end-user facility's policies and procedures concerning biohazardous materials and waste.

#### Storage

Store at controlled room temperature.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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